




# GeneProof®

NEW PRODUCT  
LAUNCH January 2020

## GeneProof MH/UU/UP Multiplex PCR Kit



### MORE INFORMATION WITH LESS SAMPLES

- PCR Kit detects and distinguishes the important STI pathogens in a single reaction tube:

*Mycoplasma hominis*

*Ureaplasma urealyticum*

*Ureaplasma parvum*

### UNIVERSAL PCR PROFILE

- Easily combinable with other PCR kits from the STI group, especially with GeneProof CT/NG/MG Multiplex PCR Kit
- Simplifies laboratory workflow

### CONTAMINATION ELIMINATION

- Ready-to-Use Master Mix contains Uracil-DNA glycosylase (UNG) and dUTPs eliminating contamination with amplification products



COMPATIBLE WITH A WIDE  
RANGE OF REAL-TIME PCR  
DEVICES



CERTIFIED  
DIAGNOSTIC TEST