



Sexually  
Transmitted  
Infections

# GeneProof®

## NEW PRODUCT LAUNCH January 2020

# GeneProof Trichomonas vaginalis PCR Kit



### HIGH SPECIFICITY

- Secured by multicopy sequence of beta-tubulin gene
- Ensures excellent sensitivity reaches up to 0.25 cp/ $\mu$ l

### RELIABLE DETECTION

- Simple sensitive detection compared to wet mount microscopy and culture

### UNIVERSAL PCR PROFILE

- Easily combinable with other PCR kits from the STI group, especially with GeneProof CT/NG/MG Multiplex PCR Kit
- Simplifies laboratory workflow

### CONTAMINATION ELIMINATION

- Ready-to-Use Master Mix contains Uracil-DNA glycosylase (UNG) and dUTPs eliminating contamination with amplification products



COMPATIBLE WITH A WIDE  
RANGE OF REAL-TIME PCR  
DEVICES



CERTIFIED  
DIAGNOSTIC TEST

[www.geneproof.com](http://www.geneproof.com)

GUARANTEED CONTROL OVER THE COMPLETE PROCESS  
OF DEVELOPMENT, MANUFACTURING AND DISTRIBUTION  
OF ALL THE OFFERED GENEPROOF PRODUCTS