



GeneProof[®]

NEW PRODUCT
LAUNCH January 2020

GeneProof Treponema pallidum PCR Kit



HIGHLY SPECIFIC DETECTION

- Secured by highly conservative polA gene for DNA polymerase I

CONTAMINATION ELIMINATION

- Ready-to-Use Master Mix contains Uracil-DNA glycosylase (UNG) and dUTPs eliminating contamination with amplification products

UNIVERSAL PCR PROFILE

- Simplifies laboratory workflow
- Easily combinable with other PCR kits from the STI group, especially with GeneProof Herpes Simplex Virus (HSV-1/2) PCR Kit
- Suitable for differential diagnostics of genital ulcer, appearing in the primary stage of syphilis



COMPATIBLE WITH A WIDE
RANGE OF REAL-TIME PCR
DEVICES



CERTIFIED
DIAGNOSTIC TEST